

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Adenosine 5'-triphosphate-Agarose - lyophilized powder

Product Number:

A2767

MDL:

MFCD00801332

Storage Temperature:

-20 °C

**TEST****Specification**

Appearance (Color)

White to Off-White

Appearance (Form)

Powder

Ligand Immobilized

1 - 5

(micromoles ATP/mL packed gel)

Swelling Volume (ml/g)

Conforms

10.0 to 15.0 ml per gram

Specification: PRD.1.ZQ5.10000009748

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.