

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

4-Aminobenzoic acid - ReagentPlus®, ≥99%

Product Number:

A9878

CAS Number:

150-13-0

MDL:

MFCD00007894

Formula:


C₇H₇NO₂

Formula Weight:

137.14 g/mol

Storage Temperature:

2 - 8 °C



TEST

Specification

Appearance (Color)

White to Off-White

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Purity (Titration by NaOH)

≥ 99 %

Purity (TLC)

≥ 99 %

Solubility (Color)

Colorless to Faint Yellow

Solubility (Turbidity)

Clear

50 mg/ml, EtOH

Recommended Retest Period

4 Years

Specification: PRD.0.ZQ5.10000044706

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.