


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Nitro-1,8-naphthalic anhydride - 95%

Product Number: 324345
CAS Number: 6642-29-1
MDL: MFCD00013446
Formula: C₁₂H₅NO₅
Formula Weight: 243.17 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared Spectrum
Purity (HPLC)

Light Yellow to Dark Yellow to Dark Beige
Powder or Crystal or Chunk(s)
Conforms to Structure
≥ 94.5 %

Specification: PRD.2.ZQ5.10000051520

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.