

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

## **Product Specification**

NaF

Product Name:

Sodium fluoride - BioXtra, ≥99%

Product Number: \$7920 CAS Number: 7681-49-4

MDL: MFCD00003524

Formula: FNa

Formula Weight: 41.99 g/mol

TEST Specification

Appearance (Color) White to Off-White Appearance (Form) Powder

Solubility (Color) Colorless
Solubility (Turbidity) Clear

1 M, H2O Insoluble matter  $$\leq$ 0.1 \%$  Chloride (Cl) Pass

Chloride (CI) P < /= 0.005%

Sulfate Pass

< 0.03%</pre>
Aluminum (Al)
< 0.001 %</pre>

 Lead (Pb)
 ≤ 0.001 %

 Iron (Fe)
 ≤ 0.003 %

 Copper (Cu)
 < 0.0005 %</td>

 Zinc (Zn)
 < 0.0005 %</td>

 Calcium (Ca)
 < 0.0005 %</td>

 Magnesium (Mg)
 < 0.0005 %</td>

 Potassium (K)
 < 0.005 %</td>

Specification: PRD.1.ZQ5.10000009427

%Purity

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.