

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Sodium hydroxide - anhydrous, free-flowing, Redi-Dri™, reagent grade, ≥98%, pellets

Product Number:

901915

CAS Number:

1310-73-2

NaOH

Formula:

NaOH

Formula Weight:

40.00 g/mol

**TEST****Specification**

Appearance (Color)

White

Appearance (Form)

Powder or Pellets

Solubility (Color)

Colorless

Solubility (Turbidity)

Clear

100 mg/ml, H<sub>2</sub>O

Titration with HCl

≥ 98 %

Impurity

≤ 1.0 %

As Na<sub>2</sub>CO<sub>3</sub>

Identity by ICP - Sodium

Positive

Recommended Retest Period

2 years

Specification: PRD.1.ZQ5.10000088757

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.