


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Aldicarb-(N-methyl-¹³C,₃d₃ carbamoyl-¹³C) - ≥98 atom %, ≥98% (CP)

Product Number: **733865**
CAS Number: 1261170-75-5



Formula: 13C₂C₅D₃H₁₁N₂O₂S
Formula Weight: 195.27 g/mol
Storage Temperature: -20 °C

TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Purity (HPLC)	≥ 98.0 %
Proton NMR Spectrum	Conforms
Carbon NMR Spectrum	Conforms
Mass Spectrometry	Conforms
Isotopic Enrichment	≥ 98 %

Specification: PRD.0.ZQ5.10000102172

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.