

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2-Methyl-3-(3,4-methylenedioxyphenyl)-propanal - ≥98%

Product Number:

W523909

CAS Number:

1205-17-0

MDL:


MFCD00067053

Formula:

C₁₁H₁₂O₃

Formula Weight:

192.21 g/mol



TEST

Specification

Appearance (Color)

Colorless to Dark Yellow

Appearance (Form)

Liquid

Refractive index at 20 °C

1.532 - 1.536

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 98.0 %

Arsenic (As)

≤ 3.0 ppm

Cadmium (Cd)

≤ 1.0 ppm

Mercury (Hg)

≤ 1.0 ppm

Lead (Pb)

≤ 10.0 ppm

Specification: PRD.2.ZQ5.10000015237

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.