

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Poly(ethylene glycol) diacrylate - average M_n 250

Product Number:


475629

CAS Number:

26570-48-9

MDL:

MFCD00081876



Storage Temperature:

2 - 8 °C

TEST	Specification
Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Liquid
Proton NMR spectra	Conforms to Structure
Molecular Number	248 - 268
Water (by Karl Fischer)	≤ 1.00 %
Acid Value	≤ 1.00
(mg KOH/g)	
Color Test	≤ 200 APHA
Inhibitor	80 - 120 ppm
(MEHQ)	
Purity	≥ 92.0 %

Specification: PRD.1.ZQ5.10000033972

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.