

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Hydroxy-TEMPO - 97%

Product Number: 176141
CAS Number: 2226-96-2

Formula: C₉H₁₈NO₂
Formula Weight: 172.24 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Colour)	Orange
Appearance (Form)	Powder or Crystals
Carbon Content	60.6 - 65.0 %
Nitrogen Content	7.8 - 8.4 %
Infrared Spectrum	Conforms to Structure

Specification: PRE.0.ZQ5.10000017822

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.