

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

4-Chlorobenzenesulfonyl chloride - 97%

Product Number:

133698

CAS Number:


98-60-2

Formula:

C₆H₄Cl₂O₂S

Formula Weight:

211,07 g/mol



TEST

Specification

Appearance (Colour)

White to Off White

Appearance (Form)

Powder or Crystals or Solid or Chunks

Purity (GC)

≥ 96,5 %

Infrared Spectrum

Conforms to Structure

Specification: PRE.0.ZQ5.10000020030

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.