

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

holo-Transferrin bovine - Iron-saturated, BioReagent, suitable for cell culture

Product Number:

T1283

CAS Number:

11096-37-0

MDL:

MFCD00131337

Storage Temperature:

2 - 8 °C

TEST**Specification**

Appearance (Color)

Orange to Brown

Appearance (Form)

Powder

Solubility (Color)

Orange

Solubility (Turbidity)

Clear

50 mg/mL, H₂O

Protein by Biuret

≥ 95 %

Iron (Fe)

≥ 0.1 %

Identity

Pass

Of Bovine Origin

Agarose Electrophoresis

97 - 100 %

Endotoxin Assay

≤ 30 EU/mg

Cell Culture Test

Pass

Recommended Retest Period

6 years

Specification: PRD.0.ZQ5.10000054207

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.