

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Ammonium acetate - reagent grade, ≥98%

Product Number:**A7262**

CAS Number:

631-61-8

MDL:

MFCD00013066

Formula:


C₂H₇NO₂

Formula Weight:

77.08 g/mol

Storage Temperature:

2 - 8 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
1 g/ml, H ₂ O	
Titration with NaOH	≥ 98 %
Recommended Retest Period	-----
2 years	

Specification: PRD.1.ZQ5.10000025105

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.