


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
2-Methyltetrahydrofuran – BioRenewable, anhydrous, ≥99%, Inhibitor-free

Product Number: **673277**  
CAS Number: 96-47-9  
MDL: MFCD00005367



Formula Weight: 86.13 g/mol

TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 99.0 %
Residue on Evaporation	≤ 0.0003 %
Peroxide (ppm)	≤ 20 ppm
AS H2O2	
Water (by Karl Fischer)	≤ 0.002 %
Note	-----

For 100 mL unit size, water content may exceed above specification; not to exceed 0.005%

Specification: PRD.2.ZQ5.10000012259

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.