


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Clonazepam – powder

Product Number: C1277
CAS Number: 1622-61-3
MDL: MFCD00057746
Formula: C₁₅H₁₀ClN₃O₃
Formula Weight: 315.71 g/mol



TEST

Specification

Appearance (Color)	White to Light Yellow
Appearance (Form)	Powder
Solubility (Color)	Yellow
Solubility (Turbidity)	Clear
20 mg/mL, Acetone	
Infrared spectrum	Conforms to Structure
Purity (TLC)	≥ 98 %
Recommended Retest Period	-----
14 years	

Specification: PRD.0.ZQ5.10000081902

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.