

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Chitosan - from shrimp shells, practical grade

Product Number:


417963

CAS Number:

9012-76-4

MDL:

MFCD00161512



TEST	Specification
Appearance (Color)	White to Beige
Appearance (Form)	Powder or Flakes
Solubility (Color)	Faint Yellow
Solubility (Turbidity)	Slightly Hazy to Hazy
10 mg/ml, 1 M Acetic acid	
Deacetylation	≥ 75 %

Specification: PRD.2.ZQ5.10000045474

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

