

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

## **Product Specification**

**Specification** 

Product Name:

**TEST** 

Glial Cell Line-derived Neurotrophic Factor human - recombinant, expressed in E. coli, lyophilized powder, suitable for cell culture, ≥98% (SDS-PAGE and HPLC)

Product Number: G1777

Storage Temperature: -20 ℃

Formulation	Conforms
Lyophilized from 0.2 um-filtered solution of 10 mM	
sodium citrate and 150 mM NaCl, with 0.5 mg BSA per ug GDNF	
Bioactivity	Conforms
ED50: 5-10 ng/mL by Dopamine uptake by rat ventral	
mesencephalic cultures	
Purity	Conforms
> or = 98% by SDS and HPLC	
Endotoxin Assav	Conforms

Specification: PRD.0.ZQ5.10000046117

< or = 0.1 ng/ug GDNF

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.