

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

Boron phosphate 99.995%

Product Number:**255955**BPO₄

CAS Number:

13308-51-5

Formula:

BO₄P

Formula Weight:

105.78 g/mol

TEST	Specification
Appearance (Color)	White
Appearance (Form) Powder and/or Chunks	Conforms to Requirements
X-Ray Diffraction	Conforms to Structure
ICP Major Analysis Confirms Boron and Phosphorus Components	Confirmed
Trace Metal Analysis	≤ 60.0 ppm
Purity 99.995% Based On Trace Metals Analysis	Meets Requirements

Remarks:

Specification Date : 12/17/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.