

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Methyl 6-methylpyridine-3-carboxylate - 97%

Product Number:

284777

CAS Number:

5470-70-2

MDL:


MFCD00006340

Formula:

C₈H₉NO₂

Formula Weight:

151.16 g/mol

**TEST****Specification**

Appearance (Color)

Conforms to Requirements

Colorless to Amber (Liquid) or
White to Tan or Orange (Solid)

Appearance (Form)

Solid or Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 96.5 %

Specification: PRD.1.ZQ5.10000013911

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.