

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

SSC Buffer 20× Concentrate – for Northern and Southern blotting, solution

Product Number:

S6639

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Solution
pH	6.9 - 7.1
at 25 degrees Celcius	
DNAse, Exonuclease Detection	None Detected
RNAse Detection	None Detected
NICKase, Endonuclease Detection	None Detected
Suitability	Suitable
Recommended Retest Period	Complete
1 year	

Specification: PRD.2.ZQ5.10000037107

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.