

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

RNase B Glycoprotein Standard from bovine pancreas - Proteomics Grade

Product Number: R1153
CAS Number: 9001-99-4
MDL: MFCD00132184

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Suitability	Suitable
Suitable as a substrate with PNGase F for N-linked deglycosylation	
Purity by SDS-PAGE	≥ 90 %

Specification: PRD.0.ZQ5.10000028840

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.