


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Poly(methyl methacrylate) – average M_w ~15,000 by GPC, powder

Product Number: 200336
CAS Number: 9011-14-7
MDL: MFCD00134349
Formula: $[\text{CH}_2\text{C}(\text{CH}_3)(\text{CO}_2\text{CH}_3)]_n$



TEST	Specification
Appearance (Color)	White
Appearance (Form) Powder or Crystalline Powder and/or Lumps	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Loss on Drying	≤ 3 %
Solubility (Turbidity) 25mg/ml in Dichloromethane	Clear to Hazy
Solubility (Color)	Colorless to Faint Yellow

Specification: PRD.0.ZQ5.10000018831

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.