

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lithium azide solution - 20 wt. % in H₂O

Product Number:

480525

MDL:

MFCD00041956

Formula:

LiN₃

Formula Weight:

48.96 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Titration by HCL	15 - 25 % wt
ICP Major Analysis	Confirmed
Confirms Lithium Component	

Specification: PRD.1.ZQ5.10000029482

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.