

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Anti-Mouse IgG (Fab specific)–FITC antibody produced in goat – affinity isolated antibody, buffered aqueous solution
Product Number: F5262

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Yellow
IgG Concentration (OD 280nm)	3.0 - 6.5 mg/ml
F/P Molar Ratio	2.5 - 6.5
IF Cells	Human periferal blood lymphocytes
IF-Titer	1:128-1:256
IHC - Tissue	Human Appendix or Human Tonsil
IHC - Fixation	FFPE-Formalin Fixed Paraffin Embedded
IHC - Titer	1:150-1:300
Recommended Retest Period	-----
3 Years	

Specification: PRE.0.ZQ5.10000016617

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.