

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Nitric Oxide Synthase, Universal antibody produced in rabbit – affinity isolated antibody, buffered aqueous solution

Product Number: **N217**

MDL: MFCD01324876

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
PBS containing 1 mg/mL BSA and 0.05% sodium azide	
Immunohistochemistry	Conforms
Working dilution 1:100	
Western Blot	Conforms
Working dilution 1:1000	
Concentration	Conforms
0.14 mg/mL	

Specification: PRD.0.ZQ5.10000059407

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.