

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Dipeptidyl Peptidase IV human – recombinant, expressed in baculovirus infected Sf9 cells, pkg of
≥1.0 units/vial, ≥10 units/mg proteinProduct Number: **D4943**

Storage Temperature: -20 °C

TEST	Specification
µg protein/ml	≥ 10
Protein by Bradford	
units/mg Protein (Bradford)	≥ 10
Unit Definition: One unit will produce 1.0 micromole of P-Nitroaniline from Gly-Pro-P-Nitroanilide per minute at pH 7.6 at 37°C.	

Specification: PRD.3.ZQ5.10000019206

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.