


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Oxalic acid - purified grade, 99.999% trace metals basis

Product Number: 658537
CAS Number: 144-62-7
MDL: MFCD00002573
Formula: C₂H₂O₄
Formula Weight: 90.03 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Purity (Titration by NaOH)	98.0 - 102.0 %
Purity	Meets Requirements
99.999% Based On Trace Metals Analysis	
Trace Metal Analysis	≤ 15.0 ppm

Specification: PRD.2.ZQ5.10000025868

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.