


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(4-Chlorophenyl)acetonitrile - 96%

Product Number: C28006
CAS Number: 140-53-4
MDL: MFCD00001918
Formula: C₈H₆ClN
Formula Weight: 151.59 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
(Run as Melt) AJP 03-20-19
Assay

White to Yellow
Solid
Conforms to Structure

≥ 95.5 %

Specification: PRD.1.ZQ5.10000024955

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.