


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Water - for HPLC

Product Number: 270733
CAS Number: 7732-18-5
MDL: MFCD00011332
Formula: H₂O
Formula Weight: 18.02 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Colorless
Appearance (Form)	Liquid
UV Absorbance 400nm	≤ 0.005
UV Absorbance 250nm	≤ 0.005
UV Absorbance 210nm	≤ 0.01
UV Absorbance 205nm	≤ 0.01
Residue on Evaporation	≤ 0.0003 %
HPLC Absorbance 254 nm	≤ 2 mAU
HPLC Absorbance 210 nm	≤ 10 mAU
Expiration Date Period	-----
12 Months	

Specification: PRD.5.ZQ5.10000003870

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.