

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tris-Borate-EDTA buffer - 5x Concentrate

Product Number: **T6400**

MDL: MFCD00236358



TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Solution
pH	8.2 - 8.4
Tested at 25°C	
Total Aerobic Microbial Count	
Report Result	
RNAse Detection	None Detected
Purity by Titration	
Report Result	
Suitability	Suitable
For use in Gel Electrophoresis after dilution to working concentration, and no DNAse/Nickase activity detected.	
Recommended Retest Period	Complete
1 year	

Specification: PRD.3.ZQ5.10000068262

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.