


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-(Dimethylamino)benzaldehyde - 98%

Product Number: 109762
CAS Number: 100-10-7
MDL: MFCD00003381
Formula: C₉H₁₁NO
Formula Weight: 149.19 g/mol



TEST

Specification

Appearance (Color)

White to Yellow and Faint Beige to Beige

Appearance (Form)

Faint Green to Green

Powder or Crystals or

Infrared spectrum

Crystalline Chunk(s) or Granule or Flakes

Purity (Titration by HClO₄)

Conforms to Structure

Purity (GC)

97.5 - 102.5 %

≥ 97.5 %

Specification: PRD.0.ZQ5.10000026352

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.