

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Activated charcoal - DARCO®, 12-20 mesh, granular

Product Number: 242241
CAS Number: 7440-44-0
MDL: MFCD00133992
Formula: C
Formula Weight: 12.01 g/mol

C

TEST	Specification
Appearance (Color)	Black
Appearance (Form) Small Chunks	Conforms to Requirements
Particle Size Distribution + 12 Mesh	≤ 10 %
Particle Size Distribution + 10 Mesh	≤ 3 %
Particle Size Distribution -20 Mesh	≤ 5 %
Moisture Content	≤ 12 %
Dust Content Dust Content	≤ 0.40 %
Vendor Information Darco Activated Carbon Is A Registered Trademark Of American Norit Company, Inc.	Confirmed

Specification: PRD.0.ZQ5.10000015213

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.