


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Nitroaniline - 98%

Product Number: N9829
CAS Number: 99-09-2
MDL: MFCD00007782
Formula: C₆H₆N₂O₂
Formula Weight: 138.12 g/mol



TEST	Specification
Appearance (Color) Yellow, Yellow-Orange or Gold	Conforms to Requirements
Appearance (Form) Powder, Crystals or Fibers and/or Chunks	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Purity (Titration by HClO ₄)	97.5 - 102.5 %
Purity (GC)	≥ 97.5 %

Remarks:

Reviewing Chemist :

Specification: PRD.0.ZQ5.10000040752

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.