

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium tetrachloropalladate(II) - 98%

Product Number:**205796**

CAS Number:

10025-98-6

MDL:


MFCD00011373

Formula:

Cl₄K₂Pd

Formula Weight:

326.43 g/mol



TEST	Specification
------	---------------

Appearance (Color)

Conforms to Requirements

Brown or Brown-Green

Appearance (Form)

Conforms to Requirements

Crystals or Crystalline Powder

X-Ray Diffraction

Conforms to Structure

Palladium (Pd)

31.8 - 33.4 %

Specification: PRD.0.ZQ5.10000036157

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.