

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sodium nitrite - ACS reagent, ≥97.0%

 Product Number:
 237213

 CAS Number:
 7632-00-0

 MDL:
 MFCD00011118

Formula: NNaO2 Formula Weight: 69.00 g/mol

NaNO₂

TEST Specification

Appearance (Color) Conforms to Requirements

White to Off-White to Pale Yellow

Appearance (Form) Powder or Crystals X-Ray Diffraction Conforms to Structure Titration by KMNO4 \geq 97.0 %

% NaNO2

Insoluble matter $\,\leq\,$ 0.01 %

C= 10%, H2O
Chloride Content < 0.005 %

Sulfate $$\leq$ 0.01 \%$$ Calcium (Ca) $$\leq$ 0.01 \%$$ Iron (Fe) $$\leq$ 0.001 \%$

Heavy Metal < 0.001 % as Lead

Meets ACS Requirements Current ACS Specification

Recommended Retest Period ------

2 Years

Potassium (K)

Specification: PRD.3.ZQ5.10000028670

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

< 0.005 %