


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tributylborane solution – 1.0 M in THF

Product Number: 179906



MDL: MFCD00009423
Formula: C₁₂H₂₇B
Formula Weight: 182.15 g/mol

TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Mass Balance	0.95 - 1.10 M

Specification: PRD.0.ZQ5.10000065064

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.