

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Lipodisq™ - Styrene:Maleic Anhydride Copolymer 3:1, Pre-hydrolyzed

Product Number: **L9170**
CAS Number: 26762-29-8

TEST	Specification
Appearance (Color)	Off White to Yellow
Appearance (Form)	Powder or Flakes
IR Spectrum	Conforms to Structure
% Composition by size exclusion chromatography(SEC)	98 - 100
Specific Gravity 1.05-1.34 at 25 Deg C	Conforms
pH 5% solution in H2O, < 9	Conforms

Specification: PRD.1.ZQ5.10000037183

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.