


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3,5-Pyrazoledicarboxylic acid monohydrate - 97%

Product Number: P56801
CAS Number: 303180-11-2
MDL: MFCD00149323
Formula: C₅H₄N₂O₄ · H₂O
Formula Weight: 174.11 g/mol



TEST

Specification

| | |
|---------------------|-------------------------------|
| Appearance (Color) | White |
| Appearance (Form) | Powder or Crystal or Chunk(s) |
| Infrared spectrum | Conforms to Structure |
| Titration with NaOH | 96.5 - 103.5 % |
| Purity (HPLC) | ≥ 96.5 % |

Specification: PRD.1.ZQ5.10000004517

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.