


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
5-Bromo-2,4-dimethoxybenzoic acid - 97%

Product Number: 593788
CAS Number: 32246-20-1
MDL: MFCD05865172
Formula: C₉H₉BrO₄
Formula Weight: 261.07 g/mol



TEST

Specification

Appearance (Color)	White to Yellow
Liquid: Colorless to Yellow	
Solid: White to Off-White	
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Purity (Titration by NaOH)	96.5 - 103.5 %
GC (area %)	≥ 96.5 %

Specification: PRD.0.ZQ5.10000037353

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.