

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lanthanum - ingot (under oil), 99.9% trace rare earth metals basis

Product Number:**261130**

CAS Number:

7439-91-0

MDL:

MFCD00011066

Formula:

La

Formula Weight:

138.91 g/mol

La

TEST	Specification
------	---------------

Appearance (Color)

Conforms to Requirements

Dark Grey

Appearance (Form)

Conforms to Requirements

Ingots in Oil

ICP Major Analysis

Confirmed

Confirms Lanthanum Component

Purity

Conforms

99.9% Based On Rare Earth Analysis

Trace Rare Earth Analysis

≤ 1500.0 ppm

Specification: PRD.0.ZQ5.10000052999

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.