

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tetracaine hydrochloride – meets USP testing specifications

Product Number: **T3937**
 CAS Number: 136-47-0
 MDL: MFCD00038912
 Formula: C₁₅H₂₄N₂O₂ · HCl
 Formula Weight: 300.82 g/mol



TEST	Specification
Identity	Pass
Water (by Karl Fischer)	≤ 2.0 %
Residue on ignition (Ash)	≤ 0.1 %
Organic Impurities	Pass
Assay	98.0 - 102.0 %
Anhydrous basis	
Residual Solvents	Meets Requirements
Expiration Dated Product	-----

Specification: PRD.2.ZQ5.10000029704

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.