


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-[(Trimethylsilyl)ethynyl]benzaldehyde - 97%

Product Number: 523380
CAS Number: 77123-57-0
MDL: MFCD02093765
Formula: C₁₂H₁₄OSi
Formula Weight: 202.32 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Dark Yellow and Faint Orange to Dark Orange and Faint Brown to Dark Brown
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Carbon	68.7 - 73.7 %
Purity (GC)	≥ 96.5 %

Specification: PRD.0.ZQ5.10000056841

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.