

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium selenocyanate – reagent grade, 97%

Product Number:**216186****KSeCN**

CAS Number:

3425-46-5

MDL:

MFCD00011366

Formula:

CKNSe

Formula Weight:

144.08 g/mol

TEST**Specification**

Appearance (Color)

Conforms to Requirements

White to Gray or Tan

Appearance (Form)

Solid

ICP Major Analysis

Confirmed

Confirms Potassium and Selenium Components

Purity (Titration by HClO₄)

≥ 52.6 %

% Selenium

Specification: PRD.0.ZQ5.10000018349

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.