


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Copper(II) sulfate – ReagentPlus®, ≥99%

Product Number: C1297
CAS Number: 7758-98-7
MDL: MFCD00010981
Formula: CuO4S
Formula Weight: 159.61 g/mol



TEST	Specification
Appearance (Color)	White to Light Grey
Appearance (Form)	Powder
Solubility (Color)	Green
Solubility (Turbidity)	Clear to Hazy
100 mg/ml, 5 M HCl	
Loss on Drying	≤ 4.5 %
Titration with Na ₂ S ₂ O ₃	≥ 99 %
Dry basis	

Specification: PRD.1.ZQ5.10000016728

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.