


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium salicylate – ReagentPlus®, ≥99.5%

Product Number: S3007
CAS Number: 54-21-7
MDL: MFCD00002440
Formula: C7H5NaO3
Formula Weight: 160.10 g/mol



TEST

Specification

Appearance (Color)	White to Off White
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Solubility (Color)	Colorless to Very Faint Yellow
Solubility (Turbidity)	Clear
50 mg/ml, H ₂ O	
Titration with HClO ₄	99.5 - 100.5 %
Recommended Retest Period	-----
4 Years	

Specification: PRD.1.ZQ5.10000035566

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.