

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

Sodium formate reagent grade, 97%

Product Number:**107603**

CAS Number:

141-53-7

MDL:


MFCD00013101

Formula:

CHNaO2

Formula Weight:

68.01 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Conforms to Requirements
Powder, Crystals or Granules or Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Purity (Titration by HCL04)	96.5- 103.5 %
Recommended Retest Period	-----
3 Years	

Remarks:

Reviewing Chemist :

Specification Date : 06/21/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.