

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
DT-2 trifluoroacetate salt - ≥95% (HPLC), film

Product Number: **D0444**

Formula: C₁₂₂H₂₂₃N₅₃O₂₃ · xC₂HF₃O₂
Formula Weight: 2,800.42 g/mol
Storage Temperature: -20 °C

TEST	Specification
Identification	Conforms
Structural	
Elemental Composition	Conforms
C ₁₂₂ H ₂₂₃ N ₅₃ O ₂₃	
Purity (HPLC)	≥ 95 %
Appearance (Form)	Conforms
lyophilized powder	

Specification: PRD.0.ZQ5.10000068120

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.