

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1,3-Bisbenzyl-2-oxoimidazolidine-4,5-dicarboxylic acid - 97%

Product Number:

650730

CAS Number:

59564-78-2

MDL:


MFCD03701609

Formula:

C₁₉H₁₈N₂O₅

Formula Weight:

354.36 g/mol



TEST

Specification

Appearance (Color)

White to Yellow

Appearance (Form)

Crystal to Powder

Proton NMR spectrum

Conforms to Structure

Purity (Titration by NaOH)

≥ 95.0 %

Specification: PRD.0.ZQ5.10000052835

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.