

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

L-Lactic Dehydrogenase from bovine muscle - Type X, ammonium sulfate suspension, ≥600 units /mg protein

Product Number: L1378
CAS Number: 9001-60-9
MDL: MFCD00131466

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

mg protein/ml (Biuret)	3.0 - 25.0
units/mg protein	≥ 600
Pyruvate Kinase	≤ 0.01 %

Specification: PRD.0.ZQ5.10000048873

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.