

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1-Butanol - for molecular biology, ≥99%

Product Number:

B7906

CAS Number:

71-36-3

MDL:

MFCD00002964

Formula:

C₄H₁₀O

Formula Weight:

74.12 g/mol

**TEST****Specification**

Appearance (Clarity)

Clear

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Proton NMR Spectrum

Conforms to Structure

Purity (GC)

≥ 99.0 %

Water (by Karl Fischer)

≤ 0.05 %

Suitability

Suitable

Suitable for use in extraction of nucleic acid solutions

Specification: PRD.1.ZQ5.10000008722

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.