


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Biotin picolyl azide

Product Number: 900912

Formula Weight: 622.74 g/mol
Storage Temperature: -20 °C

TEST	Specification
Appearance (Color)	White to Beige
Appearance (Form)	Powder or Crystals
¹ H NMR Spectrum	Conforms to Structure
Carbon Content	46.6 - 57.5 %
Nitrogen Content	16.0 - 20.0 %
Sulfur Content	4.6 - 5.7 %
Specific Rotation	10.0 - 25.0 °
0.5%, DMSO	
Purity (HPLC)	≥ 89.5 %

Specification: PRD.0.ZQ5.10000099725

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.